

EXAM TABLE PAPER ROLL

FEATURES

- Smooth texture surface table paper for greater patient comfort
- Premium table paper for higher quality and durability
- Prevent cross-contamination from patient to patient
- Safe material: 100% virgin wood pulp paper
- Cutting Line - easy to tear
- Protect exam table from dirt and moisture, helping it to stay clean and last longer
- Suitable for chiropractic exam or massage
- Work with exam table or massage table paper holder, save the space



QUALITIES







Durable, Breathable, Absorbent, Less Noise and Less Slipping, Strong, Tear Resistant, Eco-friendly, Comfortable & Soft



SPECIFICATION

Type	Single-use, disposable examination paper roll		
Material	Wood pulp paper	Sterile	No
Color	White	Packing	12 Rolls per Carton
Standard Width	50 cm	Standard length	20m, 22m & 24m (other lengths on request)
Core Diam.	50mm	Certificates	ISO13485 & ISO9001

APPLICATIONS

	Medical Consultation & Examination Room		Physiotherapy & Rehabilitation Centre
	Dental Clinics		Industrial / Factory operations & anti-oil pollution
	Spa & Aesthetic Centres		General Housekeeping,