



# Examination Paper roll

Disposable Roll is mainly used in Hospitals, Beauty Salons, Nursing Homes, Hotels, and so on to cover examination tables or massage beds.

When patients come into direct contact with medical furniture, the potential transfer of harmful contaminants poses a risk to subsequent patients who may be exposed to these pathogens. In such instances, the use of Medical Exam Table Paper becomes imperative, serving as a protective measure to guard medical furniture against bacteria, germs, and contagions. Duracare's Exam Table Paper Rolls play a crucial role in preventing cross-contamination between patients.

Designed to cover treatment chairs, examination table, massage tables, chiropractic chairs, and other medical equipment or furniture that comes into direct contact with patients, our Exam Table Paper provides an effective protective barrier. It acts as a shield between the patient and the treatment table, offering a crucial defense against the transfer of contaminants from one patient to another through contact with common surfaces. Duracare's Exam Table Paper ensures a protective barrier against harmful microorganisms, enhancing the overall safety and hygiene standards in medical settings.



- ✓ 100% Safe and Sanitary.
- ✓ Easy to use and gets rid of unnecessary laundry.
- ✓ Made with Wood pulp paper.
- ✓ Colors: White
- ✓ Protect users from bacteria and contagions
- ✓ Latex-free.

SPECIFICATION			
Type	Single-use, disposable examination paper roll		
Material	Wood pulp paper	Sterile	No
Color	White	Packing	12 Rolls per Carton
Standard Width	50 cm	Standard length	20m, 22m & 24m (other lengths on request)
Core Diam.	50mm	Certificates	ISO13485 & ISO9001



M-17, Phase 1 - Dubai Industrial City - DIC - Dubai

[www.duracaremedical.com](http://www.duracaremedical.com)

[info@duracaremedical.com](mailto:info@duracaremedical.com)

[sales@duracaremedical.com](mailto:sales@duracaremedical.com)

+971 4 323 7848

+971 56 488 1881



**DURACARE**