

Technical Data Sheet

DISPOSABLE 3 PLY FACE MASK WITH EAR LOOP

FEATURES

- Soft, comfortable, and breathable
- Large Coverage protection
- Pliable nose piece for a secure fit
- Ear loop style face mask
- 3Ply nonwoven material

Also available with head loop









QUALITY STANDARDS

Complies with:
 EN14683:2019+AC:2019 TYPE II: Type \geq 98 %: Bacterial Filtration Efficiency
 ISO 22609:2004: Test method for resistance against penetration by synthetic blood 80mmHg
 EN ISO13485:2016: Medical devices - Quality management systems - Requirements for regulatory purposes



APPLICATIONS

	Medical Purpose, Examination & Laboratory applications		Industrial Purpose & PPE
	Healthcare & Nursing		General Housekeeping
	Food & Beverage Industry & Hospitality		IT Industry



الإمارات
THE EMIRATES

Made in UAE