

Technical Data Sheet

DISPOSABLE UNDERPADS


FEATURES

- Multi-Purpose Protection
- Clean, hygienic & Reliable Incontinence Solution
- Effective Hygiene Control
- Quickly Absorbs Moisture
- Waterproof penetration
- Skin friendly & Comfortable Fluff
- Latex Free



PRODUCT STRUCTURE





- 1st layer: Soft non-woven film
 2nd layer: Tissue paper
 3rd layer: Fluff pulp mixed with SAP, to absorb the liquids
 4th layer: Tissue paper
 5th layer: PE film, to prevent any leakage through



SPECIFICATION

Type	Disposable Underpad – Nursing Pads		
Material	PE, Wood Pulp & SAP	Colors	Blue - White - Green
Standard Sizes	600x600 - 600x900 (Custom size on demand)	Certificates	ISO13485 & ISO 9001
Packing	25 pcs per bag / 150pcs per carton	Sterile	No

APPLICATIONS

	Hospital & Medical facilities applications		Pet training
	Healthcare & Nursing		Baby Care



الامارات
THE EMIRATES
Made in UAE